

NEGATIVE PRESSURE PUMP

# NewEraI<sup>®</sup>



## NewEraI<sup>®</sup> CANISTER

MODEL NUMBER: 16670-NEWERA



- Compact
- Foldout bed rail hanger
- Intuitive touchscreen user interface
- Carrying bag included
- IV pole adaptor available
- Reliable pump motor





The NewEra-I is a compact negative pressure pump that is controlled through an intuitive touchscreen interface for ease of operation. Alarms for pressure tolerance and high vacuum as well as to alert for low battery, tube blockage, Inactive pump, and when the canister is full help personnel stay attuned to the patient's needs. A carry bag and foldout bed rail hangers help during patient transport.



### NewEra I PUMP SPECIFICATIONS

REQUIREMENT	SPECIFICATION	NOTES
Portable	Yes	Carry bag available
Weight	2.43 lbs.	
Pressure Range	-10 TO -200 mmHg	Adjustable in 5 digit increments
Battery Life		Up to 10 hours
Cycles	Yes	Intermittent or Continuous
Displays	Digital	Displays pressure setting, mode of operation and alarm activated
	Status LED	Operational LED
	Other LED	Slow Mode controls rate of negative pressure increase or decrease
Auto Lockout	Yes	User controlled time settings
Blue Tooth Data Capture	Yes	Alarms, On/Off, runtimes, notifications
Pressure Tolerance Alarm	Yes	If the vacuum is from 5 to 80 mmHg below the user controlled set point
High Vacuum Alarm	Yes	If the vacuum is from 5 to 80 mmHg higher than the user controlled set point
Other Alarms	Yes	Low battery, Tube blockage, Inactive pump, canister full
Alarm Mute Button	Yes	
Main Power Option	Yes	
Power Input	Power Adapter 5V - 2.0Amp	

### NewEra I CANISTER SPECIFICATIONS

REQUIREMENT	SPECIFICATION	NOTES
Canister Sizes	900ml	
Solidifier with Canister	Yes	
Fluid Measurement Marks	Yes	
Single touch canister removal	Yes	
Antibacterial Filter	Yes	3 micron bacterial filter with a 32 mm diameter to prevent microbial contamination
Tubing Clamps	Yes	

Integrated Healing Technologies™